

AED Checklist

AED location		Check date	
AED make		AED serial no.	
Checked by		Signature	

Item	Checked	Comments/action required	
AED Location			
Visible and unobstructed			
AED sign displayed			
AED Condition			
Check battery status/in date			
Check service indicator light			
AED exterior undamaged			
AED Components			
Two sets of AED pads			
Two sets of paediatric pads*			
Pads in date and in sealed pack			
Check expiry date on pads			
Pocket mask present			
Examination gloves present			
Safety Razors present			
Towel or material to dry patient			
Clothing shears			
Faults or defects noted			
Reported to:		Date:	

** If applicable*